

SHIPPING KITS

When traveling parents need to send breast milk home, our shipping kits offer a reliable solution. We always overnight ship—because milk matters.

PUMP & SHIP

This is our bestselling shipping solution. It's designed for pumping parents who need to ship refrigerated breast milk while away from home.

✓ Two size options

- 34 oz capacity
 - 72 oz capacity
- Includes 6 storage bags
- Includes 12 storage bags

✓ 60 hours of cooling

Keeps breast milk safely refrigerated from pickup to delivery.

✓ Includes everything

Pre-labeled FedEx Priority Overnight return label, reusable cooler, tamper-proof seals, and full instructions.

✓ Door-to-door pickup

Parents can schedule a FedEx pickup from their hotel, office, or conference center..



Pump & Ship takes the guesswork out of breast milk logistics and gives parents peace of mind while they're away. It's a favorite for business trips, trainings, and travel-heavy roles.

FREEZE & SHIP

Built for frozen milk. Whether a parent is on a multi-day trip or just prefers to freeze their milk before sending it home, the Freeze & Ship is a reliable, parent-friendly option.

✓ Two size options

- 48 oz capacity
- 90 oz capacity

✓ 96 hours of cooling

Maintains frozen temperatures throughout the shipping process.

✓ Includes everything

Each kit comes with a ready-to-use cooler, tamper-evident seals, detailed packing instructions, and a pre-addressed FedEx Priority Overnight shipping label.

✓ Flexible delivery options

Parents can schedule a pickup or drop it at the nearest FedEx location—whatever works best for their day.



The Freeze & Ship is especially helpful for those on longer trips.



CARRY-ON COOLERS

Ideal for travel days, conferences, and shorter trips, our carry-on coolers provide a convenient solution for transporting breast milk.

We provide TSA guidance with every cooler to assist parents during security screenings.

PUMP & TOTE

Designed for parents who prefer to carry their milk with them, the Pump & Tote offers:

- ✓ **Two size options**
 - 34 oz capacity  Includes 6 storage bags
 - 72 oz capacity  Includes 12 storage bags
- ✓ **60 hours of cooling**

Ensuring breast milk remains fresh during extended trips.
- ✓ **TSA-compliant design**

Facilitates smooth passage through airport security.
- ✓ **Easy activation**

Features a push-button cooler for immediate use.



Each kit comes with a tote bag for discreet and convenient transport.

FREEZE & TOTE

For parents needing to transport frozen milk, the Freeze & Tote provides:

- ✓ **96 hours of cooling**

Maintains breast milk in a frozen state throughout travel.
- ✓ **Compact and TSA-friendly**

Designed for hassle-free carry-on.
- ✓ **48 oz capacity**

Accommodates various breast milk storage bag types.



This cooler is perfect for long flights or weekend trips.

CHECK-IN COOLERS

Suitable for longer journeys or when carrying larger quantities of breast milk, our check-in coolers are designed for checked luggage.



Includes clear
TSA guidance for
a hassle-free
experience.

PUMP & CHECK / FREEZE & CHECK

Designed for parents who prefer to carry their milk with them, the Pump & Tote offers:

- ✓ **108 oz capacity**
Typically holds a 3-5 day supply of breast milk.
- ✓ **Minimum of 90 hours refrigeration**
Keeps breast milk at the optimal temperature during transit.
- ✓ **Comprehensive kit**
Includes storage bags, circular seals, tamper-proof stickers, travel bag, and luggage tag.

This solution is ideal for extended trips where carrying breast milk on board isn't feasible.



HOW IT WORKS

Ordering and using Milk Stork products is a straightforward process:



Employee Access

Log in using the company-provided credentials.



Product Selection

Choose the appropriate solution based on travel plans and breast milk transportation needs.



Delivery Options

Have the kit shipped to a preferred location—home, office, hotel, or other destinations



Usage

Follow the provided instructions to pack and transport breast milk.



Support

Access customer service for assistance throughout the process.

No codes. No hoops. Just support when it's needed most.