

How it Works



Pump according to your regular nursing schedule and refrigerate your milk



Pack. Activate the cooling unit in your Milk Stork cooler (instructions included). Then, pack your refrigerated breast milk



Ship. Schedule a FedEx pickup on your own, or ask your hotel to help. Your breast milk will be overnighted home*



Tote. Carry your cooler of milk home with you in the reusable tote provided**

Disclaimers:

* FedEx Priority Overnight's shipping capabilities allow for next-day delivery to most destinations. FedEx does not operate on Sundays. Alaska, Hawaii and rural destinations may have longer transit times. Occasionally, weather and other events can cause disruptions to FedEx's service. Milk Stork cannot control for these disruptions. Milk Stork's Pump & Ship coolers provide at least 72hr of refrigeration.

** Whether carried on or checked as luggage, your Milk Stork cooler AND your breast milk will be inspected by TSA and may also be X-rayed. If you're toting your breast milk internationally, it is the user's responsibility to research and comply with each country's, airport's, and security agency's regulations for traveling with breast milk. Milk Stork is not responsible for inconsistent or incorrect enforcement of breast milk travel regulations and/or inspection of the Milk Stork coolers in the United States or internationally.

Customer Service Information

✉ | info@milkstork.com

☎ | 888-207-6909

🕒 | M-Sa, 6am - 8pm CT



milkstork[®]

Travel solution for breastfeeding moms

Milk Stork is the first and only breast milk delivery service for business traveling moms, providing no-fuss, refrigerated, express shipping or easy toting of breast milk home to baby.

www.milkstork.com/stryker

Express love[®]

The Pump & Ship

For a mom who needs a no-fuss solution to overnight her refrigerated breast milk home.

- Pharmaceutical-grade shipping cooler with easy-to-use push button activation
- Provides at least 72 hours of refrigeration
- Shipped to mom's travel destination
- Pre-labeled and pre-addressed with FedEx Priority Overnight shipping labels for next day delivery
- Also includes: Breast milk storage bags, packing instructions and shipping seals

Available in Two Sizes:

34 oz: For a mom who would like to overnight a daily supply of breast milk.

72 oz: For a mom who would like to overnight more milk less often.

The Pump & Tote

For a mom who would prefer to carry her refrigerated breast milk home.

- Pharmaceutical-grade cooler with easy-to-use push button activation
- The size of an airplane carry-on bag
- Provides up to 60 hours of refrigeration
- Also includes: Breast milk storage bags, packing and toting instructions, and a tote bag
- Shipped to mom's destination and ready to tote home

Available in Two Sizes:

34 oz: For a mom who would like to carry approximately a day's supply of breast milk home with her. Weighs 6.3 lbs when full.

72 oz*: For a mom who would like to carry a multi-day supply of breast milk home with her. Weighs 9.6 lbs when full.

*72 oz. is over two days' worth of milk for most moms.

