

We Set the Gold Standard for Breast Milk Shipping

Milk Stork is the first breast milk delivery service for business traveling moms, providing no-fuss express shipping or easy toting of refrigerated breast milk home to baby.

How it Works



Pump

According to your regular nursing schedule and refrigerate your milk.



Pack

Activate the cooling unit in your Milk Stork cooler. Then, pack your refrigerated milk in the breast milk storage bags provided.



Ship

All shipments come pre-labeled with a FedEx Priority Overnight label. Schedule a FedEx pickup on your own or ask your hotel to help.



Tote / Check

Carry your cooler of milk home, or if traveling internationally, check it as luggage. Use the Milk Stork travel bags provided.



For more information
or to place your order, go to:
portal.milkstork.com/stryker

DOMESTIC

The Pump & Ship

Available in three Sizes:

34oz

72oz

- 12" X 5" X 8"
- 72 hr minimum of refrigeration
- One day's worth of milk for most moms.

- 12" X 8" X 8"
- 72 hr minimum of refrigeration
- Over two days' worth of milk for most moms.

108oz

- 12" X 8" X 10"
- 72 hr minimum of refrigeration
- 3/4 days' worth of milk for most moms.

The Pump & Tote

Available in two Sizes.

34oz

72oz

- Weighs 6.3 lbs when full.
- 12" X 8" X 8"
- 60 hr maximum of refrigeration.
- One day's worth of milk for most moms.

- Weighs 9.6 lbs when full.
- 12" X 8" X 8"
- 60 hr maximum of refrigeration.
- Over two days' worth of milk for most moms.

Disclaimers:

*FedEx Priority Overnight's shipping capabilities allow for next-day delivery to most destinations. FedEx does not operate on Sundays. Alaska, Hawaii and rural destinations may have longer transit times. Occasionally, weather and other events can cause disruptions to FedEx's service. Milk Stork cannot control for these disruptions. Milk Stork's Pump & Ship coolers provide at least 72hr of refrigeration.

**Whether carried on or checked as luggage, your Milk Stork cooler AND your breast milk will be inspected by TSA and may also be X-rayed. If you're toting your breast milk internationally, it is the user's responsibility to research and comply with each country's, airport's, and security agency's regulations for traveling with breast milk. Milk Stork is not responsible for inconsistent or incorrect enforcement of breast milk travel regulations and/or inspection of the Milk Stork coolers in the United States or internationally.



Questions?

support.milkstork.com | +1-510-356-0221 Mon-Sat, 6am - 8pm CT

DOMESTIC

The Freeze & Tote

Available in one Size:

48 oz

- 12" X 8' X 8"
- Weighs 5lbs when full
- Keeps breast milk frozen for a minimum of 96 hours without dry ice
- One day's worth of milk for most moms.

The Freeze & Ship

Available in two Sizes:

48 oz.

90oz

- 12" X 8' X 8"
- Weighs 5 lbs when full
- Keeps breast milk frozen for a minimum of 96 hours without dry ice
- One day's worth of milk for most moms.

- 12" X 8' X 10"
- Weighs 9.5 lbs when full
- Keeps breast milk frozen for a minimum of 96 hours without dry ice
- 3/4 day's worth of milk for most moms.

The Freeze & Check

Available in one Size:

90oz

- 12" X 8' X 10"
- Weighs 9.5 lbs when full
- Keeps breast milk frozen for a minimum of 96 hours without dry ice
- 3/4 day's worth of milk for most moms.



DOMESTIC

The Cool Commuter

- External: 10.5 x 8 x 7 inches
- Internal: 8.5 x 5.25 x 5 inches
- Holds up to 24 oz. of breast milk and keeps it refrigerated for 15 hours
- Includes 2 ice packs



INTERNATIONAL

The International Pump & Check

Check your refrigerated breast milk as luggage.

Available in One Size:

108oz

- Weighs 16 lbs when full.
- 12" X 8" X 10"
- 90 hr minimum of refrigeration.
- Over three days' worth of milk for most moms.

Disclaimers:

*FedEx Priority Overnight's shipping capabilities allow for next-day delivery to most destinations. FedEx does not operate on Sundays. Alaska, Hawaii and rural destinations may have longer transit times. Occasionally, weather and other events can cause disruptions to FedEx's service. Milk Stork cannot control for these disruptions. Milk Stork's Pump & Ship coolers provide at least 72hr of refrigeration.

**Whether carried on or checked as luggage, your Milk Stork cooler AND your breast milk will be inspected by TSA and may also be X-rayed. If you're toting your breast milk internationally, it is the user's responsibility to research and comply with each country's, airport's, and security agency's regulations for traveling with breast milk. Milk Stork is not responsible for inconsistent or incorrect enforcement of breast milk travel regulations and/or inspection of the Milk Stork coolers in the United States or internationally.

Questions?

support.milkstork.com | +1-510-356-0221 Mon-Sat, 6am - 8pm CT