



Moms on a Mission

portal.milkstork.com/stryker

How it Works



According to your regular nursing schedule and refrigerate your milk.



Activate the cooling unit in your Milk Stork cooler. Then, pack your refrigerated milk in the breast milk storage bags provided.



All shipments come pre-labeled with a FedEx Priority Overnight label. Schedule a FedEx pickup on your own or ask your hotel to help.

OR



Carry your cooler of milk home, or if traveling internationally, check it as luggage. Use the Milk Stork travel bags provided.

Milk Stork is the first breast milk delivery service for business traveling moms, providing no-fuss, refrigerated, express shipping or easy toting of breast milk home to baby.

#MomOnAMission